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EN

SUMMARY AND EXPLANATION

Diagnostic kit ST 16 is standardized identification system for common species Gram positive catalase positive microorganisms, especially of *Staphylococci*. Set is based on 16-17 miniaturized biochemical tests and software database. List of all microorganisms, for which is the kit determined, is placed at the end of leaflet.

PRINCIPLE

Kit consists of 16 wells strip of microtitration plate in classic 96 wells format containing dehydrated substrates, where ST 16 sp is in the form of strippable microtitration plates and ST 16 fp is in the form of undivided full plate. Reconstitution of substrates runs by inoculation of bacterial suspension. During incubation occurs colour change in wells because of metabolic activity of microorganisms. Evaluation of results of tests is visual on base of colour scheme or by colour formulation described in leaflet. Results of identification are obtained from evaluation table or by help of evaluation software, which can be found at www.diagnostics.sk/idmicro.

KIT CONTAINS- ACCORDING TO PACKAGING 60 tests (sp) / 150 TESTS (fp)

- 10 / 25 microtitration plates of ST 16
- 60 /150 result forms
- 10 / 25 incubation packets
- 1 information leaflet

REQUIRED MATERIALS

Reagents:

- Unbuffered saline 3,5-5 ml
- Paraffin oil (Ref. 3001)
- NIT reagent (Ref. 3005)
- Zn (Ref. 5001)
- Optional:
 - VP and VP reagent (Ref. 2004 and 3004)
 - PYR and PYR reagent (Ref. 2003 and 3003)
 - Identification software

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Materials:

- Pipettes, tampons, loops, burner, tubes and other basic laboratory equipment.

WARNINGS AND SAFETY PRECAUTIONS

- **For in vitro diagnostics use and microbial control**
- **For professional use only**
- Follow the instructions exactly!
- Used strips should be considered as potentially infectious and this must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packaging is intact. Do not use damaged kit.

By interpretation of results anamnesis of patient, source of sample, morphology of colony, microscopic morphology of batch and if necessary, results of all overrun tests, mainly results of antibiogram have to be considered.

STORAGE CONDITIONS

Diagnostic kits are delivered in multilayer packets on base of aluminium and organic polymers. Part of each packet is silica gel desiccant. Store kits at temperatures from +2 to +25°C. Expiration date is placed on each packaging.

Put unused opened microtitration strip in packed Al packet with original silica gel desiccant, close packet carefully and store it at laboratory temperature. Opened product can be stored two weeks or to date of expiration, in case when expiration date passes earlier.

SAMPLES

Isolate microorganisms, which have to be identified, from suitable unselective cultivation medium (for example blood agar, etc.) according to standard microbiological techniques. Make Gram colouring and microscope pure culture. Make catalase test. Confirmed isolates identify on ST 16.

RECOMMENDED PROCEDURE

Preparation of inoculum

- Open the tube with 2,5-5ml of unbuffered sterile saline.
- Take same well isolated colonies by inoculation loop from 18-24 hours old culture.
- Turbidity of homogenized suspension must be 3 McF. This suspension must be used immediately after preparation!

TIP: Eventually make proof of purity by braying of same loop or tampon, by which was suspension done.

This Petri dish can be used next day to make other additive tests!

Preparation of microtitration plate

- Prepare microtitration plate
- Mark strips with numbers of examining cultures.

Inoculation

- Inoculate by 0,1 ml of well homogenized suspension into each well of strip.
- Cover tests URE to ORN (wells H1 to F1) with 2 – 3 drops of paraffin oil.

Incubation

- Put the microtitration plate to packed PE packet, then bend end of packet under plate – this avoids dehumidifying of bacterial suspension.
- Incubate by common atmosphere at temperature 35 ± 2 °C for 18 - 24 hours.

TIP: For optimal incubation increase humidity in incubator for example by adding beak of water or incubate by incubator with managed humidity.

EVALUATION AND INTERPRETATION

We highly recommend marking the result of NOV and COA! Evaluate the tests after incubation with help of evaluating form, colour scheme, or by results of control batch.

Test bGA/NIT is bifunctional so after evaluating of primary reaction another result can be obtained from evaluated well of microtitration plate.

Make following optional bifunctional test

Well A1: bGA / NIT - after evaluation of primary reaction add 1-2 drops of NIT reagent and evaluate test NIT. Confirm negative test results with Zn dust. Color change after application of Zn dust means negative result. No color change confirm positive result!

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- Write down result of bifunctional tests to form for evaluating results or to evaluating software.

IDENTIFICATION

Result of identification can be obtained by:

• Identification by identification table:

Compare results of tests (page. 5) and make identification by results of tests in identification table (page 6).

• Identification by identification software:

Enter results of individual tests. When some of the test cannot be evaluated, it can be excluded. Software allows inputting additional tests and so it increases identification effectiveness. microID software is for free in web based and web independent version.

QUALITY CONTROL

Quality of diagnostic kits is systematically controlled. Chemicals are bought only from certified companies and quality of these chemicals is confirmed by analytical certificate. The functionality of the kits is tested by collection of control strains, controlled and tested is also present of bacterial contamination. Kits are exposed to tests of higher temperatures and samples of each batch are saved for right advisement of later reclamations. For the need of own proof function, use recommend bacterial strains (page 5).

CONSTRAINTS OF METHOD AND MOST OFTEN CAUSES OF WRONG IDENTIFICATION

- Diagnostic kit ST 16 is determined for identification of bacteria's named in this leaflet only.
- Only pure culture of microorganism can be used.
- Tests were not covered by paraffin oil.
- Contamination of wells by inoculum of next strip.
- Used culture is atypical batch.
- Some point of leaflet was not kept

DETERMINATION CHARACTERISTICS

Tested were 62 collection strains and strains with clinical background, but veterinary strains mentioned in database were also tested:

- 91,9 % strains were identified correctly.
- 8,1 % strains were not identified

WASTE LIQUIDATION

Work with material as with potentially infectious agents. Liquidate rubbish according to internal procedures and according to legislative of your country. Kit components do not contain dangerous chemicals.

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